

MTXQCvX Part4: PROJECT NAME *

Secret *Secret*
Christin *Kempa Lab*

Contents

Metmax-parser for MTXQC	1
Project settings	1

MTXQC part 4 provides the transformation of Metmax-derived files for the usage as input files for MTXQC part 1. This report provides three modes - (1) basic - creates tables required for GC-MS performance, (2) Quant - Generation of ManualQuantTable (3) Inc - Calculation of ¹³C-isotope incorporation.

Keywords: MTXQCvX, pSIRM time course, cell extracts, manual validation, quantities, stable isotope incorporation

Metmax-parser for MTXQC

Project settings

```
#set path for figure export and size

set_input = "input/"
set_output = "output/"
## subfolder for postprocessing

#directory definition and figure_name definition
if (params$spath == "") {
  path_setup = ""
  set_fig = paste0(path_setup, 'figure/MTXQCp4-')
} else {
  path_setup = paste0(params$spath, "/")
  set_fig = paste0(path_setup, 'figure/MTXQCp4-')
}

knitr::opts_chunk$set(fig.width = 8, fig.align = 'center', fig.height = 7,
  fig.path = set_fig,
  echo = FALSE, #TRUE - show R code
  warning = FALSE, #show warnings
  message = TRUE,
  eval = TRUE
) #show messages

## Correct input format of files defined! metmax

## File imported! file annotation

## File imported! sample_extracts
```

*Template MTXQCvX part 4 written by Christin Zasada, Kempa Lab

```
## PeakArea matrix imported!  
## No file for m/z 73 values defined!  
## No file with MIDs defined!  
## Please define your internal standard in the conversion_metabolite.csv file!  
## Alkane standard annotation detected in conversion_metabolite.csv.
```